

International Workshop on Future Linear Colliders 2014 (LCWS14)

# Status of Beam Loss Monitoring system for the CLIC Two-Beam Module and Damping Rings

<u>M. Kastriotou</u>, M.J. Boland, E. Branger, A. Degiovanni, E. Effinger, L. Devlin, S. Doebert, W. Farabolini, D. Gutkov, E. B. Holzer, J. van Hoorne, P.D. Jackson, R.L. Lillestol, S. Mallows, J. L. Navarro, E. Nebot del Busto, R.P. Rasool, J. Schmidt, S. Sousa, F. Tecker, S. Vinogradov, C. P. Welsch, B. Wooley, and M. Zingl





- Compact Linear Collider (CLIC): Large nominal power (Drive Beam 70 MW, Main Beam 14 MW)
- $\bullet$  Loss of a small fraction of the beam  $\rightarrow$  severe impact on accelerator components

#### Beam Loss Monitors:

#### → Machine Protection

 ${\ }$  Detect potentially dangerous instabilities  $\rightarrow$  destruction of the surrounding materials, activation of the accelerator components

Prevent subsequent injections into the Main Beam Linac and the Drive Beam decelerators

#### $\rightarrow$ Beam Diagnostics

 The fraction of beam particles lost has to be carefully controlled to achieve good beam transmission

Localise and characterise the Beam Loss distribution

## Outlook

- BLM system for the Two-Beam Module (TBM):
  - Requirements
  - Detector Technology choice: Ionisation Chambers (IC), Diamonds and Quartz Crystals
  - Optical fibre based BLM (OBLM) system: Challenges
    - Radiation Hardness
    - Choice of photon sensors
  - BLMs at the Test Beam Line (TBL): Measurements and Monte Carlo model
  - Installation at the dogleg experiment and first measurements
  - Installation at the TBM and BLM program
- First look at Damping Rings
  - Considerations for detector technology choice
  - Measurements at Australian Synchrotron
  - Sensitivity studies
  - OBLM single shot calibration
  - User fill measurements (Topup injections)
  - Dynamics (CLIC-like fill scans)

#### • Summary and conclusions



## Requirements: A recall

#### • Sensitivity

Losses should be controlled to 0.1% to avoid luminosity losses due to beam loading variations

#### Dynamic range

- Detect onset (10%) of dangerous losses
- 0.1% of DB
- 0.001% of MB
- FLUKA simulations

	Sensitivity [Gy]	Dynamic Range
DB 2.4 GeV	10-7	5×10 <sup>5</sup>
DB 0.24 GeV	5×10 <sup>-8</sup>	2×10 <sup>5</sup>
MB 9 GeV (Type 1)	10-9	<b>10</b> <sup>6</sup>
MB 1.5 TeV (Type 4)	10 <sup>-8</sup>	10 <sup>5</sup>

#### A MULTI-TEV LINEAR COLLIDER BASED ON CLIC TECHNOLOGY

CLIC CONCEPTUAL DESIGN REPORT



## **Detector technology**

• A system based on Ionization Chambers satisfies the requirements of Machine Protection -(CLIC Conceptual Design Report)

- LHC (1.5 l N<sub>2</sub>): Sensitivity 5.4 x10<sup>-9</sup> C/Gy (1 pA). DR > 10<sup>+7</sup>
- But requires a large number of BLMs ( > 40000)
- Optical fibre based BLMs (OBLMs)

Pure Silica core multimode fibres as <u>Cherenkov radiators</u>. Both ends of fiber couples to photon-sensors

Potential of covering full beam line in a cost-effective way. Position resolution within 30 cm of beam loss has been demonstrated



- Fast BLMs (diagnostics oriented): Diamond and Quartz:
  - Diamonds (1 ns): LHC for bunch by bunch beam loss monitoring
  - Quarzt crystals + photon sensor: time response depends on volume of crystal
  - (~50ps/cm<sup>3</sup>) and photon sensor. Potential for ultra-fast Beam Loss Monitoring

# Cherenkov fibers: Challenges I

#### Radiation hardness

- Material, manufacturer, type of radiation
- Pure Si core with high OH recommended
- Radiation Induced Attenuation strongly dependent on  $\lambda$
- SiO<sub>2</sub> fibers rather insensitive for 800 nm and above.



#### J. Kuhnhenn

Fiber irradiation test at Fraunhofer institute 10kGy@0.22Gy/s

# Cherenkov fibers: Challenges II

- Choosing a photon sensor
  - Quantum Efficiency
  - Saturation Effects
  - Time Response
  - Achievable Bandwidth





Avalance Photo Diode (APD) Photon Multiplying tube (PMT) Multi Pixel Photon Counter (MPPC)

#### maria.kastriotou@cern.ch

#### LCWS14, 7th October 2014

## Cherenkov fibers: Challenges II



# BLMs at the Test Beam Line: Cherenkov fibre

- Pure Silica 200µm core multimode optical fibre located 28cm on top of the beam line
- 25 (75) meters of read-out fibre in the downstream (upstream) direction
- 14400 pixel (25um) 3x3mm<sup>2</sup> Hamamatsu MPPCs
- Trans-impedance amplifier (opAmp THS3601 GBW 300 MHz) with Rf = 500  $\Omega$







#### Hamamatsu MPPC

# Localised BLMs at TBL

- Four detector technologies
  - Ionisation chamber x2
  - Aluminium Cathode EMission (ACEM) + PMT x8
  - pCVD diamond x1

PEP-II (Cherenkov crystal + Fast PMT) x 1

Read via 250 MHz 8 Bit ADC (OASIS)
 Occasionally 5GHz/10Gs/s for timing investigations





Installation

during 2013

# Localised BLMs at TBL

- Four detector technologies
  - Ionisation chamber x 2
  - Little Ionisation Chamber (LIC) x 4
  - pCVD diamond
  - PEP-II (Cherenkov crystal + Fast PMT ) x 1
  - ACEM (Al foil + PMT) x 1
- Read via 250 MHz 8 Bit ADC (OASIS)
  Occasionally 5GHz/10Gs/s for timing investigations





Installation

during 2014

#### Monte Carlo Model losses at the TBL



#### Beam loss measurements at the TBL



	Dotoctor	<b>FLUKA Simulation (pC/e)</b>		Measurement	<b>Combined Beam</b>
	Detector	Systematic variation	Most realistic value	(nC)	(mA)
fiber -	Downstream fiber	(2.5 - 19.0) x 10 <sup>-6</sup>	9.0 x 10 <sup>-6</sup>	2.36	260
	Upstream fiber	(0.037 - 4.4) x 10 <sup>-6</sup>	1.4 x 10 <sup>-6</sup>	0.15	370
Locali- sed BLM	ACEM	(0.93 - 5.0) x 10 <sup>-20</sup>	2.3 x 10 <sup>-20</sup>	1.21	34
	PEP-II	(1.1 - 11.0) x 10 <sup>-20</sup>	4.2 x 10 <sup>-20</sup>	1.25	13
	pCVD	(0.97 - 4.0) x 10 <sup>-20</sup>	2.1 x 10 <sup>-20</sup>	0.42	19

#### ✓ The BLMs detect losses lower than the BPM sensitivity!

maria.kastriotou@cern.ch

## Installation at the dogleg experiment

- Measurements of the RF Breakdown and dark current in the accelerating structure
- Comparison of the signals without and with beam loading
- Localised
  - pCVD diamond
    40 dB Broadband
    Current Amplifier
- OBLM
  - Ø 900 µm fibre +
    SiPM
  - Custom-made lowpass voltage filter and crate modules



# Preliminary data without beam

- System optimised for dark current measurements
- RF pulse length and shape
- → Under improvement for RF Breakdown diagnostics.





# Installation at the Two-Beam Module (TBM)

#### • 4 Little Ionisation Chambers

- Downstream of quads at the Drive Beam
- Downstream of accelerating structures at the Main Beam
- <u>OBLMs</u> along the Drive Beam and the Main Beam
- Supports for the integration on the TBM.

# <image>

#### TBM BLM aim of studies

- Nominal beam losses at the TBM
- Comparison with simulations
- Extrapolations to high energy
- X-talk of losses





# **CLIC damping rings**

#### Considerations

- Synchrotron light production as a means of beam cooling  $\Rightarrow$  insensitivity of BLM system
- Sensitive to potential "quench inducing" beam losses

#### • OBLM system

- The full ring can be covered with ~4 BLMs
- Cherenkov-based fibre BLMs are insensitive to synchrotron radiation

#### •Tested in the Australian Synchrotron (AS)

- LINAC (10 m): 90 keV to 100 MeV
- Booster (130 m): 100MeV to 3 GeV
- Storage Ring (216 m): 3 GeV (Similar parameters and flexibility)



Parameter	AS	CLIC DR
Energy (GeV)	3.0	2.86
Intensity (elec)	$9.0 \cdot 10^{+11}$	$1.28 \cdot 10^{+12}$
Number of bunches	300	312
Pulse lenght (ns)	600	156
Circ. lenght (m)	216	427.5
$f_{rev}$ (MHz)	1.38	0.73
Bunch spacing (ns)	2	0.5
$\gamma \epsilon_x$ (nm rad)	58708	472
$\gamma \epsilon_u$ (nm rad)	< 5	4.8

## Australian Synchrotron measurements

- Installation of several prototypes in the AS
- Two (one) 7 m (5 m) optical fibres with 365  $\mu m$  (200  $\mu m$ ) SiO2 core:
  - Multi Pixel Photon Counter (MPPC)
  - Photon Multiplier Tube (PMT)
  - Avalanche Photon Diode (APD)
- Pin diode, Nal and NE102 scintillators in neighboring locations for comparison









## **OBLM** calibration

- Intensity variation
  10<sup>3</sup>-10<sup>9</sup> charges
- Single bunch directed onto closed scraper



- OBLM system sensitive to localized losses down to ~10<sup>-9</sup> and linear up to 0.001 of the full CLIC pulse.
  - Saturation effect observed for MPPC due to small spot illumination

# BLM measurements during topup injections

- AS works on topup mode to keep a constant 200 mA for > 24 hours
- Frequency Analysis
  - SR ~ 0.5 mA / Booster ~ 1.1mA / Injection efficiency ~80%
  - harmonic number Booster/SR = 216/360
  - 1.25x10<sup>+9</sup> electrons lost within the first several turns
- Frequency Analysis
  - 11 kHz band (synchrotron tune)
  - 1.38 MHz (rev frequency)
  - 400 kHz band (vertical tune)







# **CLIC-like** beam

• CLIC-like bunch length, 150 ns

•11 independent skew quadrupole settings to change the coupling and vertical emittance

8 pulses per setting

 Measurement of beam losses during betatron coupling scan





Amount of losses seems to increase steadily with the emittance

11

- Different detector technologies examined for the CLIC Two-Beam Module BLMs
- Optical fibre BLM development: ongoing tests at CTF3
- The ideal BLM system will combine OBLMs and localised detectors for beam diagnostics
- OBLM system tested on damping ring like conditions
  - ✓ Measure Beam losses 10<sup>3</sup>- 10<sup>9</sup> electrons
  - Indications of insensitivity to SR
  - Measurement of losses for a CLIC bunch train

## Thank you!!!

# Cherenkov fibers: light spectrum

• For a quantitative estimation of beam losses we need to consider:

Cherenkov light
 spectrum (λ<sup>-2</sup>)
 Length and λ

dependence of light attenuation

 Quantum efficiency of selected photosensor (~400-900 nm)



# **OBLM** calibration

#### • Intensity variation

• Vgun increase to reduce the bunch charge down to  $\sim 10^{+5}$  charges

• Closings slits in Booster to SR to further reduce down to  $\sim 10^{+3}$  charges





# Detector technology

#### • A system based on Ionization Chambers satisfies the requirements of Machine

Protection – (CLIC Conceptual Design Report)

- LHC (1.5 | N<sub>2</sub>): Sensitivity 10<sup>-9</sup> Gy (1 pA). DR > 10<sup>+7</sup>
- LIC (0.5 | N<sub>2</sub>): Sensitivity 10<sup>-8</sup> Gy (4 pA)
- LHC-like Current to Frequency Converter: DR=10<sup>+5</sup>
- But requires a <u>large number of BLMs</u> ( > 40000)

 Electron (ion) total collection/rise time is ~500 ns/10 ns (400 μs/10 μs). No time information of origin of beam loss

#### • Optical fibre based BLMs (OBLMs)

Pure Silica core multimode fibres as <u>Cherenkov radiators</u>. Both ends of fiber couples to photon-sensors

 Potential of covering full beam line in a cost-effective way. Better time resolution (position resolution within 30 cm of beam loss has been demonstrated)



- Fast BLMs (diagnostics oriented): Diamond and Quartz:
  - Diamonds (1 ns): LHC for bunch by bunch beam loss monitoring
  - Quarzt crystals + photon sensor: time response depends on volume of crystal (~50ps/cm<sup>3</sup>) and photo sensor. Potential for ultra-fast Beam Loss Monitoring