

ATF2 IPBPM meeting

jitRun33 resolution study

N. Blaskovic, T. Bromwich, R. Ramjiawan, S. Yang

Outline

Analysing with the different calibration constants and the same sample ranges and triggers for consistency check with Siwon's reported results.

Jitter run: jitRun33_0dB_0.95_ipbpm_160316 (triggers 1:200)

Calibrations: IPAyCal1_0dB_0.95_ipbpm_160315
IPByCal1_0dB_0.95_ipbpm_160315
IPCyCal1_0dB_0.95_ipbpm_160315
IPAxCal1_20dB_1_ipbpm_160315
IPBxCal1_20dB_1_ipbpm_160315
IPCxCal1_20dB_1_ipbpm_160315

Sample numbers: 53:63 (integration)

Reference sample: 50

Resolution using IPA

Analysis results with geometric
and multi-parameter fits.
Coefficients for all parameters.

Fitting:

- A: YI and constant
- B: YI, YQ and constant
- C: YI, YQ, Ref Y and constant
- D: YI, YQ, XI, XQ and constant
- E: YI, YQ, XI, XQ, Ref Y, Ref X and constant
- F: YI, YQ, XI, XQ, Ref Y, Ref X, self XI, XQ and constant

Siwon's measured resolution
with 13 parameter fit
= 11.2634 nm

Parameter	Geometric	A	B	C	D	E	F
No. param	2	3	5	6	9	11	13
Y1I	1.464	1.068	1.225	3.395	1.386	3.431	3.393
Y2I	-0.464	-0.700	-0.359	0.372	-0.566	0.306	0.390
Const		3.422	-3.002	-60.093	-3.480	-59.830	-60.009
Y1Q			0.015	1.761	0.101	1.767	1.723
Y2Q			-0.741	-0.836	-0.659	-0.823	-0.896
YREF				-0.004		-0.077	0.299
X1I					-0.139	-0.048	-0.080
X2I					-0.054	0.328	0.411
X1Q					1.239	0.241	0.254
X2Q					0.016	-0.004	-0.004
XREF						0.000	0.000
X3I							-0.273
X3Q							0.338
Resolution (nm)	89.3	49.2	32.3	13.9	30.3	13.1	12.6
Res_s (nm)	2.5	2.7	1.8	0.8	1.7	0.7	0.7

Resolution using IPB

Analysis results with geometric
and multi-parameter fits.
Coefficients for all parameters.

Fitting:

- A: YI and constant
- B: YI, YQ and constant
- C: YI, YQ, Ref Y and constant
- D: YI, YQ, XI, XQ and constant
- E: YI, YQ, XI, XQ, Ref Y, Ref X and constant
- F: YI, YQ, XI, XQ, Ref Y, Ref X, self XI, XQ and constant

Siwon's measured resolution
with 13 parameter fit
= 11.6025 nm

Parameter	Geometric	A	B	C	D	E	F
No. param	2	3	5	6	9	11	13
Y1I	0.683	0.643	0.661	0.655	0.652	0.648	0.647
Y2I	0.317	0.633	0.156	0.082	0.162	0.087	0.105
Const		-0.973	0.455	2.307	0.440	2.279	1.821
Y1Q			-0.651	-0.621	-0.646	-0.616	-0.617
Y2Q			0.291	0.344	0.280	0.341	0.328
YREF				0.000		-0.021	0.039
X1I					-0.033	0.022	0.026
X2I					0.019	-0.004	-0.002
X1Q					0.159	-0.163	-0.140
X2Q					-0.147	0.000	0.000
XREF						0.000	0.000
X3I							-0.090
X3Q							-0.130
Resolution (nm)	89.3	80.6	11.6	10.6	10.6	10.1	10.0
Res_s (nm)	3.6	4.3	0.6	0.6	0.6	0.5	0.5

Resolution using IPC

Analysis results with geometric
and multi-parameter fits.
Coefficients for all parameters.

Fitting:

- A: YI and constant
- B: YI, YQ and constant
- C: YI, YQ, Ref Y and constant
- D: YI, YQ, XI, XQ and constant
- E: YI, YQ, XI, XQ, Ref Y, Ref X and constant
- F: YI, YQ, XI, XQ, Ref Y, Ref X, self XI, XQ and constant

Siwon's measured resolution
with 13 parameter fit
= 12.6368 nm

Parameter	Geometric	A	B	C	D	E	F
No. param	2	3	5	6	9	11	13
Y1I	-2.156	-0.529	-0.913	-0.848	-1.092	-0.976	-0.975
Y2I	3.156	0.795	2.030	-0.856	2.321	0.217	0.216
Const		3.597	-1.918	65.660	-3.413	44.097	44.089
Y1Q			0.769	1.719	0.929	1.511	1.511
Y2Q			0.596	-1.959	0.582	-1.162	-1.163
YREF				0.005		0.105	0.113
X1I					0.087	-0.032	-0.050
X2I					-0.130	-0.711	-0.711
X1Q					-1.028	0.956	0.949
X2Q					2.346	0.004	0.004
XREF						-0.001	-0.001
X3I							-0.006
X3Q							0.058
Resolution (nm)	89.3	40.5	26.5	15.2	15.8	11.5	11.5
Res_s (nm)	1.1	2.2	1.4	0.8	0.9	0.6	0.6

Conclusions




- Resolution results for IPA, IPB and IPC using a 13-parameter fit are close to Siwon's results. Differences are likely due to different fitting methods in determining calibration constants, and/or different flyer data cuts.
- Fit parameters seem very different from Siwon's fit parameters. Need to look at this in more detail to understand why.
- Also need to look at stability of parameter coefficients across different jitter runs.

Appendix

- I and Q waveforms for IPA, IPB and IPC from jitRun33

On the following slides:

- Calibration runs and fit plots for the calibration constants



jitRun33_0dB_0.95_ipbpm_160316

Number of triggers: 200
Number of samples: 190
IPBPM analysed: AY
Waveform starts at sample: 54
Reference maximum at sample: 50
IQ maximum at sample: 60


jitRun33_0dB_0.95_ipbpm_160316

Number of triggers: 200
Number of samples: 190
IPBPM analysed: BY
Waveform starts at sample: 52
Reference maximum at sample: 50
IQ maximum at sample: 59

jitRun33_0dB_0.95_ipbpm_160316

Number of triggers: 200
Number of samples: 190
IPBPM analysed: CY
Waveform starts at sample: 53
Reference maximum at sample: 50
IQ maximum at sample: 60

IPBPM Calibration



IPAyCal1_0dB_0.95_ipbpm_160315

Number of triggers: 551
 Number of samples: 190
 IPBPM analysed: AY
 Waveform starts at sample: 50
 Reference maximum at sample: 49
 IQ maximum at sample: 58


CHARGE
 $0.9418 \pm 0.0005 \times 10^{10}$

MOVERS
 Number of steps: 11
 Triggers per step:
 50 50 50 50 50 50 50 50 50 50 51
 Steps used: 1:11

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 2
 3-sigma ICT cuts: 0
 Total triggers cut: 11
 1 6 7 10 22 23 34 40
 40 47 537 538

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : -0.128 ± 0.006
 Scale factor k : 1.72 ± 0.03

IPBPM Calibration



IPByCal1_0dB_0.95_ipbpm_160315

Number of triggers: 551
 Number of samples: 190
 IPBPM analysed: BY
 Waveform starts at sample: 50
 Reference maximum at sample: 50
 IQ maximum at sample: 58


CHARGE
 $0.941 \pm 0.0005 \times 10^{10}$

MOVERS
 Number of steps: 11
 Triggers per step:
 50 50 50 50 50 50 50 50 50 50 51
 Steps used: 1: 11

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 4
 3-sigma ICT cuts: 1
 Total triggers cut: 5
 99 253 273 297 298

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : -0.166 ± 0.006
 Scale factor k : 1.19 ± 0.01

IPBPM Calibration



IPCyCal1_0dB_0.95_ipbpm_160315

Number of triggers: 551
 Number of samples: 190
 IPBPM analysed: CY
 Waveform starts at sample: 51
 Reference maximum at sample: 49
 IQ maximum at sample: 59


CHARGE
 $0.9421 \pm 0.0005 \times 10^{10}$

MOVERS
 Number of steps: 11
 Triggers per step:
 50 50 50 50 50 50 50 50 50 50 51
 Steps used: 1: 11

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 5
 3-sigma ICT cuts: 2
 Total triggers cut: 7
 66 133 163 173 188 189 231

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : 0.47 ± 0.01
 Scale factor κ : 0.84 ± 0.01

IPBPM Calibration



IPAxCal1_20dB_1_ipbpm_160315

Number of triggers: 121
 Number of samples: 190
 IPBPM analysed: AX
 Waveform starts at sample: 52
 Reference maximum at sample: 50
 IQ maximum at sample: 57


CHARGE
 $0.926 \pm 0.001 \times 10^{10}$

MOVERS
 Number of steps: 4
 Triggers per step:
 30 30 30 31
 Steps used: 1: 4

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 1
 3-sigma ICT cuts: 0
 Total triggers cut: 1
 30

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : 0.025 ± 0.006
 Scale factor k : -0.0427 ± 0.0004

IPBPM Calibration



IPBxCal1_20dB_1_ipbpm_160315

Number of triggers: 121
 Number of samples: 190
 IPBPM analysed: BX
 Waveform starts at sample: 54
 Reference maximum at sample: 50
 IQ maximum at sample: 58


CHARGE
 $0.9462 \pm 0.0008 \times 10^{10}$

MOVERS
 Number of steps: 4
 Triggers per step:
 30 30 30 31
 Steps used: 1: 4

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 0
 3-sigma ICT cuts: 0
 Total triggers cut: 0

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : 0.075 ± 0.006
 Scale factor κ : -0.0357 ± 0.0004

IPBPM Calibration



IPCxCa1_20dB_1_ipbpm_160315

Number of triggers: 121
 Number of samples: 190
 IPBPM analysed: CX
 Waveform starts at sample: 53
 Reference maximum at sample: 50
 IQ maximum at sample: 58

CHARGE
 $0.9485 \pm 0.0007 \times 10^{10}$

MOVERS
 Number of steps: 4
 Triggers per step:
 30 30 30 31
 Steps used: 1: 4

CUTS
 Saturated samples: 6
 Reference threshold cuts: 0
 3-sigma reference cuts: 2
 3-sigma ICT cuts: 0
 Total triggers cut: 2
 5 6

CALIBRATION
 Samples 53 to 63
 Integration
 Reference sample number: 50
 IQ rotation angle θ : 0.15 ± 0.01
 Scale factor κ : -0.036 ± 0.003